


Pose Correction System for Physical Therapy and Rehabilitation Using Computer Vision



Jithin Krishnan (A0249481W)

Master of Technology in Artificial Intelligence Systems NUS-ISS, National University of Singapore



Motivation: Why home Physical Therapy still needs a coach

Widespread form errors

- 55% of patients perform ≥ 1 drill incorrectly ¹
- Errors usually go unnoticed at home

Real clinical impact

- Mis-execution delays healing
- Re-injury risk increases when cues are missed

Most tools fall short

- Binary pass / fail feedback ⇒ no guidance
- Cloud / GPU needed ⇒ cost, latency, privacy

Lightweight way forward

- Real-time, joint-level, spoken advice on any laptop
- Summarizations and Analysis





Improving Rehabilitation:

Al-Powered Physical Therapy

Technical problem

• Laptop-only CV system, must spot 14 joint-angle errors and speak corrective cues in < 0.5s at $\ge 90\%$ accuracy.

Business upside

- ullet Cuts therapist load and follow-up visits ullet lower costs
- Scales telerehab to thousands without cloud / GPU fees

Impact for Singapore

- Shrinks wait-times for an ageing, physio-heavy population
- Backs Healthier SG & Smart-Nation tele-medicine goals







Key Contributions: RehabApp - Al Physiotherapist

Lightweight Architecture

3-head CNN + Bi-LSTM, only 3.4 M parameters

Granular Feedback

Joint-level advice on 14 angles – not just pass/fail

Built-in Safety

"Wrong-exercise" guard blocks unsafe cues automatically

Real-time & Open

handles 30 video fps on a mid-range CPU; Opensource React + FastAPI



Dataset: REHAB24-6²



Key facts	Details
Purpose	Benchmark pose-estimation & exercise-quality algorithms under real rehab conditions
Subjects & drills	10 adults (6♂ / 4♀, 25–50 y) perform 6 physiotherapy exercises Arm-abduction, Arm-VW, Push-ups, Leg-abduction, Lunge, Squats
Volume	65 recordings = 184,825 frames
Ground-truth	41-marker 3-D mocap Derived 26-joint skeleton in 3-D & 2-D Synced RGB videos
Why we chose it	Contains both visual & Motion-Capture Ground Truth (mocap GT) Balanced correct/wrong reps Multiple views → tests occlusion robustness

Exercise ID	Reps	Correct	Wrong	Frames
1	178	90	88	27 442
2	208	94	114	33 641
3	107	52	55	12 054
4	210	120	90	18 329
5	174	78	96	17 608
6	195	134	61	19 373
Total	1 072	568	504	128 447

Summary Table

video_id	repetition_nurexercise_id	person_id	first_frame	last_frame	cam17_orient mod	cap_erron exercise_subt lights_on	extra_	person_extra	person_corre	ctness
PM_001	1	1	1 130	324	front	0 right arm	0	1	0	1
PM_001	2	1	1 326	537	front	0 right arm	0	1	0	1
PM_001	3	1	1 538	731	front	0 right arm	0	1	0	1
PM_001	4	1	1 732	919	front	0 right arm	0	1	0	1
PM_001	5	1	1 920	1090	front	0 right arm	0	1	0	1
PM_001	6	1	1 1690	1940	front	0 right arm	0	1	0	0
PM_001	7	1	1 2060	2346	front	0 right arm	0	1	0	0
PM_001	8	1	1 2347	2554	front	0 right arm	0	1	0	0
PM_001	9	1	1 2558	2746	front	0 right arm	0	1	0	0
PM_001	10	1	1 2747	2961	front	0 right arm	0	1	0	0
PM_002	1	1	1 152	363	half-profile	0 right arm	0	1	0	1
PM_002	2	1	1 364	573	half-profile	0 right arm	0	1	0	1
PM_002	3	1	1 574	766	half-profile	0 right arm	0	1	0	1
PM_002	4	1	1 767	910	half-profile	0 right arm	0	1	0	1
PM_002	5	1	1 1000	1225	half-profile	0 right arm	0	2	0	0
PM_002	6	1	1 1226	1407	half-profile	0 right arm	0	2	0	0
PM_002	7	1	1 1408	1605	half-profile	0 right arm	0	1	0	0
PM_002	8	1	1 1606	1800	half-profile	0 right arm	0	1	0	0
PM_002	9	1	1 1801	1967	half-profile	0 right arm	0	2	0	0

Author-supplied frame-level segmentation





Data:


Augmentations and Pipeline

Temporal & Label Augmentations


- Sliding window T = 16 frames, stride = 8 → ×8 more samples
- Half-profile relabel
- Oblique views auto-marked incorrect (for 5 / 6 drills)

Feature Engineering

- 99-D vector / frame (33 keypoints × 3 coords)
- One-hot exercise ID (6 classes) appended once per window
- 14 joint-angle errors kept as labels for the regression head







System Architecture:

End-to-End Overview

Front-end (React, runs in browser)

- Captures 640 × 480 webcam video and runs MediaPipe Pose
 @ 30 Hz
- Builds 16 × 99 tensors + one-hot exercise ID; streams them via WebSocket
- Renders colour-coded skeleton, speaks cues (Speech-Synthesis API) and shows live KPI dashboard

Back-end (FastAPI, on same laptop)

- Hosts PoseQualityNet-KP → predicts quality, 14 joint-angle errors & exercise ID in ≈ 2 ms
- Majority-vote buffer (5 windows) suppresses flicker; "wrongexercise" guard blocks unsafe advice
- Sends compact JSON to frontend, logs KPIs & survey answers in SQLite (REST endpoint)




Model Architecture:

PoseQualityNet-KP

Block	Shape / Notes
Key-point encoder	1-D CNN 99 → 128 → 512
Temporal encoder	2-layer Bi-LSTM (256 h / dir)
Exercise embedding	6→ 64 MLP
Prediction heads	quality (2-way) • angle errors (14-D) • exercise (6-way)
Footprint	3.41 M params ≈ 13 MB





Features: Live Feedback & Analytics




Visual feedback

- Landmarks color change— sky-blue by default; top-3 joints with $|\Delta\theta| \ge 8^{\circ}$ turn amber.
- Textual feedback display green = good, red = bad form, amber = wrong exercise.
- suggestion (angle to to adjust) display in amber color



Real-time dashboard

- Updates every window (≈ 0.5 s):
- Pie chart correct / wrong
- Numeric tiles mean angle error,
 % correct
- 14-bin histogram of joint errors with the three worst joints called out.



Spoken feedback


 Browser Speech-Synthesis speaks the same message; joint-level advice (e.g. "straighten left knee and right elbow") has higher priority than generic banners.



Session summary

 Same charts from the real-time dashboard after the exercise end for the user to review and improve










Inference: Feedback Logic

Scenario	Example sentences	How they are produced
Correct form	Feedback : "You're on the right track!" Suggestion : —	 Quality head ⇒ p(wrong) < 0.60 and exercise ID agrees with the drop-down → state = good. All joint errors < 8° → corrections list is empty.
Bad form	Feedback : "You're doing it wrongly!" Suggestion (example): "Adjust your left hip, spine and right ankle properly."	 p(wrong) ≥ 0.60 and exercise matches with selected → state = bad_form Keep joints with Δθ ≥ 8°, sort by magnitude, take top 3 → e.g. {left-hip, spine, right-ankle}.
Wrong exercise	Feedback : "Wrong exercise! Looks like Push-ups." Suggestion : —	 1 Exercise head ≠ selected and conf ≥ 0.65 → state = wrong_ex 2 Skip joint-error logic entirely.





Results: Quantitative

Window-level performance (4112 test windows)

- Repetition quality: Acc 91.5 %, F1_w 0.915
- Exercise ID: Acc 99.5 %, F1_w 0.995
- Joint-angle error MAE: 4.73° across 14 DoF

Confusion-matrix takeaways


- Rep-quality: TP: 1094, TN: 879, FP: 116, FN: 68
- Exercise ID: only 5 mis-labels
- largest confusion Arm-VW ↔ Arm-Abduction
- 93 % of windows give clinically actionable (< 5°) angle errors


Efficiency & footprint

- 7500 windows $s^{-1} \approx 30$ fps end-to-end on mid-range CPU
- 3.41M parameters (< 8MB), no GPU required
- Latency: ≈ 2.5 ms (browser
 → server → browser)

Efficiency & footprint

- Mobile-ready: real-time, < 10
 MB, runs on-device
- Rich feedback: joint-level errors + wrong-exercise guard
- State-of-the-art accuracy on six diverse rehab drills





Head	Acc	F1_w	MAE (°)
Rep-quality	0.915	0.915	
Exercise ID	0.995	0.995	_
Angle error	_	_	4.73

Results: Qualitative



shoulder properly"

Excercise	Correct Form	Incorrect Form	Feedback/Suggestion
Ex1: Arm-abduction	Select Exercise: Arm abdustion State Arg error 7,6° Sab Arg error 7,6° Sab You're on the right track!	Select Exercise: Arm abduction Step Any arraw 9,1* You're doing it wrongly! Suggestion: Adjust your right shoulder, right wrist and right elbow property.	Correct form: "You're on the right track" Incorrect form: "You're doing it wrongly!" "Suggestion: Adjust your right shoulder, right wrist and right elbow properly"
Ex2: Arm-VW	Select Exercise: Arm VVV Stop Ang energy 23.77 75% You're on the right track!	Select Exercise: Ann VVV Ang erreat 13.0* Ang erreat 34% You're doing it wrongly! Suggestion: Adjust your left shoulder, right shoulder and right ankle properly.	Correct form: "You're on the right track" Incorrect form: "You're doing it wrongly!" "Suggestion: Adjust your left shoulder, right shoulder and right angle properly"
Ex3: Push-ups	Select Exercise: Push-ups Any arrac 14.11 You're on the right track!	Select Exercise: Plush-ups Stop Ang error 14.4" N correct 31% You're doing it wrongly! Suggestion: Adjust your left elbow, left shoulder and right shoulder properly.	Correct form: "You're on the right track" Incorrect form: "You're doing it wrongly!" "Suggestion: Adjust your left elbow, left shoulder and right

Results: Qualitative



properly"

Excercise	Correct Form	Incorrect Form	Feedback/Suggestion
Ex4: Leg-abduction	Select Exercise: Ling Mitouchian Ang error Ang error 15.0° 44% You're on the right track!	Select Exercise: Log abduction Step Aug error 13.2* You're doing it wrongly! Suggestion: Adjust your right elbow, right ankle and left ankle properly.	Correct form: "You're on the right track" Incorrect form: "You're doing it wrongly!" "Suggestion: Adjust your right elbow, right angle and left angle properly"
Ex5: Lunge	Select Exercise: Log lungs: Aug error 33.6° 37% You're on the right track!	Select Exercise: Ling larges Aury senser: 14.2* 34% You're doing it wrongly! Suggestion: Adjust your head, left shoulder and right ankle properly.	Correct form: "You're on the right track" Incorrect form: "You're doing it wrongly!" "Suggestion: Adjust your head, left shoulder and right angle properly"
Ex6: Squats	RehabApp: Al Physiotherapist Select Exercise: Separts Anguerner 20.8* 89% You're on the right track!	Select Exercise: Squals: Ang error: 38.8* Ang error: 84% You're doing it wrongly! Suggestion: Adjust your right wrist, right knee and right efbow properly.	Correct form: "You're on the right track" Incorrect form: "You're doing it wrongly!" "Suggestion: Adjust your right wrist, right knee and right elbow


Results: Qualitative



Safety demo – "Wrong exercise" guard

Exercise Selected: Ex1: arm-abduction "

Exercise Performed Ex2: arm-VW"



Feedback/Suggestion

"Wrong exercise! Looks like arm-abduction "



Session summary page

The pie chart, numeric tiles, and joint-error histogram from the session history, and an advice to guide the user



Advice

"You should focus on correcting your right shoulder, left shoulder, right elbow more"



Ablation Study:

Architecture Variants

Variant	Added Block(s)	Rep- F1	Ex-F1	MAE (°)	ΔRep- F1	Δ ΜΑΕ	FPS	Params
A	Baseline 1-D CNN only	0.797	0.993	8.49	-	-	9 160	0.25 M
В	+ Exercise embed (6 → 64 MLP)	0.819	0.994	6.18	↑ +2.2 pp	↓ -27 %	9 126	0.25 M
С	+ Bi-LSTM (2 × 256)	0.853	0.996	6.12	↑ +5.6 pp	↓ -28 %	7 256	3.40 M
FULL	B + C (embed + Bi- LSTM)	0.903	0.997	3.86	↑ +10.6 pp	↓ -55 %	7 383	3.41 M

Which Blocks Really Pay Off?

• Exercise context is cheap and valuable

o 0% extra params → +2pp Rep-F1 & - 2.3° MAE.

Temporal reasoning matters

 Bi-LSTM alone adds +5.6 pp Rep-F1 despite a 13× size jump.

Synergy, not mere addition

 Combining both cuts angular error in half (8.49 → 3.86°) and yields the largest quality gain, yet model still fits in ~ 11 MB.


Real-time preserved

 Even the full model sustains ≈30 FPS on CPU (7 k+ on GPU), meeting the livefeedback requirement.





System: Demo






Limitations

Self-occlusion

When a limb is hidden (e.g. arms resting on thighs, crossed legs) MediaPipe drops landmarks, triggering spurious "bad-form" flags.

Camera-pose sensitivity

Accuracy degrades if the webcam is pitched or rolled by > ±20 ° or if lighting is uneven.




Scope constraints

Model was trained on 10 healthy young adults; it is untested on elderly, post-operative or high-BMI users (risk of domain shift).

Low-resolution extremities

Single-view, single-person, six drills, no on-device personalisation or multi-person support (yet).



Future Works

Robustness

- Add a lightweight self-supervised pre-text task to denoise landmarks under occlusion / harsh lighting.
- Fuse two camera views when available to recover hidden joints and improve 3-D angle accuracy.

Personalized Coaching

- 30-second calibration routine to learn each user's neutral joint baselines.
- Dynamic tolerance tightening as the patient improves, plus on-device few-shot fine-tuning.

Ultra-light Edge Deployment

- Quantise the model to 4-bit QAT ONNX and export as a cross-platform mobile SDK (Android / iOS / WebAssembly).
- Optimise CPU scheduling to keep real-time feedback while cutting battery drain.

Clinical Validation

- Run a 6-week field study with postoperative patients (n≈30) vs. standard care; measure adherence, recovery time, and user satisfaction.
- Collect qualitative therapist feedback to refine cue phrasing.

Feature Expansion

- Double the exercise library (e.g. balance drills, gait training) and introduce multi-person support.
- Add multi-language text-to-speech and haptic cues for broader accessibility.





References

[1] Faber, M., Andersen, M. H., Sevel, C., Thorborg, K., Bandholm, T., & Rathleff, M. The majority are not performing home-exercises correctly two weeks after their initial instruction—an assessor-blinded study. PeerJ, 3:e1102, 2015. DOI:10.7717/peerj.1102.

[2] A. Černek, J. Sedmidubsky, and P. Budíková, "REHAB24-6: Physical Therapy Dataset for Analyzing Pose Estimation Methods," in *Proc. 17th Int. Conf. on Similarity Search and Applications* (SISAP), LNCS 14512, pp. 18–33, Oct. 2024, doi: 10.1007/978-3-031-75823-2_2.







THANK

YOU

